



Glenammer

Laboratory Test Sieves

# SQ Series

## Digital Shaker



## Features

SQ digital has a digital display panel with an installed digital timer, users can set up the desired operation sieving time between 0 sec to 99 mins 59 secs

The system features a pause/start option, which enables users to operate the shaker with intervals.

### SQ Series

#### Capacity

8 X 200mm Dia. (100mm, 150mm, 3" and 8") sieves plus lid and receiver.

6 X 300mm Dia. (12" 315mm ) sieves plus lid and receiver. Extended rods available.

#### Max sample weight

3Kg

#### Speed

3,000/min at 50 Hz

#### Height

140mm

#### Net W x H x D

240mm x 140mm x 240mm

#### Net Weight

20Kg

#### Sieve Motion

Electromagnetic

#### Electrical Supply

230 Volt, 1 ph, 50 Hz, Input power 0.045 kw, current 0.20 amps

#### Timer

Digital: 0sec - 99mins 59secs

#### Type

8 x 200mm dia  
(100mm, 150mm, 8")  
sieves plus lid and receiver