



Glenammer

Laboratory Test Sieves

Gem Series

Basic Shaker



Features

GEM Basic has an advanced power control system, which adjusts the amplitude of the sieving motion from 10 steps.

Users can adjust the power setting accordingly during the test or adjust the power for testing different types of materials.

GEM Basic is capable of holding test sieves that are 100mm diameter, 150mm diameter, 200mm diameter and 8" diameter.

GEM Basic has an accurate digital timing system.



Capacity

8 x 200mm Diam (100mm, 150mm, 3" and 8") sieves plus lid and receiver.
Extended rods available.

Net Weight

47Kg

Max Sample Weight

3Kg

Sieve Motion

3D Electromagnetic

Amplitude

0-3mm (Max amplitude depends on loading) Set digitally 0-9 steps.

Electrical Supply

230 Volts, 50Hz. Other voltages and frequencies are available.

Net W x H x D

400mm x 200mm x 400mm

Timer

Digital: 0sec - 99mins 59secs